
SORIN

XTRA

GB - ENGLISH	XTRA AUTOTRANSFUSION DISPOSABLE ACCESSORIES - INSTRUCTIONS FOR USE.....	3
IT - ITALIANO	XTRA - ACCESSORI MONOUSO PER AUTOTRASFUSIONE - ISTRUZIONI PER L'USO	6
FR - FRANÇAIS	ACCESSOIRES À USAGE UNIQUE POUR AUTOTRANSFUSION XTRA - INSTRUCTIONS D'UTILISATION	9
DE - DEUTSCH	XTRA EINWEGZUBEHÖR FÜR DIE AUTOTRANSFUSION - GEBRAUCHSANWEISUNG.....	12
ES - ESPAÑOL	ACCESORIOS DESECHABLES DE AUTOTRANSFUSIÓN XTRA - INSTRUCCIONES DE USO	15
PT - PORTUGUÊS	ACESSÓRIOS DESCARTÁVEIS DE AUTOTRANSFUSÃO XTRA - INSTRUÇÕES DE USO.....	18
NL - NEDERLANDS	XTRA AUTOTRANSFUSIEACCESSORIES VOOR EENMALIG GEBRUIK - GEBRUIKSINSTRUCTIES	21
SE - SVENSKA	ENGÅNGSTILLBEHÖR TILL XTRA AUTOTRANSFUSION - BRUKSANVISNING	24
DK - DANSK	ENGANGSTILBEHØR TIL XTRA AUTOTRANSFUSION - BRUGSANVISNING.....	27
NO - NORSK	ENGANGSUTSTYR TIL XTRA AUTOTRANSFUSJON - BRUKSANVISNING.....	30
FI - SUOMI	XTRA-AUTOTRANSFUUSION KERTAKÄYTTÖISET LISÄVARUSTEET - KÄYTTÖOHJEET	33
EL - ΕΛΛΗΝΙΚΑ	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΑΝΑΛΩΣΙΜΩΝ ΠΑΡΕΛΚΟΜΕΝΩΝ - ΣΥΣΤΗΜΑΤΟΣ ΑΥΤΟΜΕΤΑΓΓΙΣΗΣ XTRA	36
PL - POLSKI	JEDNORAZOWE AKCESORIA DO AUTOTRANSFUZJI XTRA - INSTRUKCJA UŻYTKOWANIA.....	39
CZ - ČESKY	PŘÍSLUŠENSTVÍ JEDNORÁZOVÉ SOUPRAVY PRO AUTOTRANSFUZI XTRA - NÁVOD K POUŽITÍ	42
РУ - РУССКИЙ	ОДНОРАЗОВЫЕ ПРИНАДЛЕЖНОСТИ СИСТЕМЫ АУТОГЕМОТРАНСФУЗИИ XTRA - ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ ..	45

CODE N.	PRODUCT DESIGNATION	INTENDED USE
04268	(BRB1) Blood Reinfusion Bag 1 liter	The device is intended to store the processed blood or blood components recovered from autotransfusion systems for subsequent reinfusion to the patient.
04269	WASTE BAG X	The device is intended to collect the fluids and cellular debris removed from the blood being processed by autotransfusion systems.
04270	Reservoir Outlet Clamp XRES B	The Reservoir Outlet Clamp XRES B is used for the setup of the Blood Collection Reservoir
04271	Reservoir outlet "Y" XRES B	The Reservoir Outlet "Y" XRES B allows for an additional access site to the Procedure Set

DESCRIPTION

The devices are single use, non toxic, non pyrogenic, supplied STERILE and packaged singly. Sterilized by ethylene oxide. The level of ethylene oxide residuals in the devices is within the limits established by national regulations in the country of use.

The accessories must be used in conjunction with the medical devices listed in the paragraph MEDICAL DEVICES TO BE USED WITH THE PRODUCT.

CONTRAINDICATIONS

There are no known contraindications when the device is used as intended.

SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

CAUTION

Indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Manufactured by

Sterile - Ethylene oxide sterilized

Non Pyrogenic

Contains PHTHALATE

Latex free

CAUTION: Federal law (USA) restrict this device to sale by or on the order of a physician

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number

Attention, see instruction for use

Attention, see instruction for use

This way up

Fragile; handle with care

Keep away from heat

Keep dry

ea unit

GENERAL SAFETY INFORMATION

The following is general safety information with the aim of advising the operator in preparing to use the device. Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

WARNING

- The device must be used in accordance with the instructions for use provided in this manual.
- The device is intended to be used by professionally trained personnel only.
- Do not use the device if it has been cracked, dropped or otherwise physically damaged.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Avoid any condition which may cause the blood temperature to exceed 37°C (98.6°F).
- Keep dry. Store at room temperature.
- Always apply and maintain a correct dose and accurate monitoring of the anticoagulant.
- For single use and for single patient use only: during use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body and due to its specific design it cannot be fully cleaned and disinfected at the end of use. Therefore the reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if STERILE.
- The device contains phthalates. Considering the nature of body contact, the limited contact duration and the number of treatments per each patient; the amount of phthalates which might be released from the device do not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorized local representative.

CAUTION

- For a detailed description of the circuits, refer to SORIN XTRA operator's manual.
- Refer to the SORIN XTRA operator's manual for complete instructions on use of the device
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician

SETUP

⚠ WARNING

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.
- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc. to prevent damage to the device.
- Make the connections using a suitable aseptic technique.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate parts of the device. This could cause damage which may compromise the integrity and proper functioning of the device.
- Due to the possibility of operator exposure to blood borne pathogens (such as HIV, hepatitis viruses, bacteria, Cytomegalovirus, etc.) when handling extracorporeal blood circuits, adequate precautions should be taken at all times to prevent the exposure to and transmission of such agents.
- Carefully observe this set for leaks before and during use. Leakage may result in loss of sterility or loss of blood and/or fluid. If leakage is observed before or during use, replace the leaking component or appropriately retighten the leaking connection.

GENERAL SETUP INFORMATION

1. Remove the product from the storage carton. Using aseptic technique, remove the external wrapping and proceed with following steps.
2. Discard after use.

04268 (BRB1) BLOOD REINFUSION BAG 1 LITER

SET UP OF THE (BRB 1) BLOOD REINFUSION BAG 1 LITER

1. Make sure the clamp on the bowl set reinfusion line is closed. Hang the BRB 1 to the upper stage of the IV pole on the autotransfusion machine.
2. Connect the female luer-lock connector on the bag line to the male luer-lock on the bowl set reinfusion line. Ensure the connection is tight.
3. Open the clamp on the bowl set reinfusion line. Make sure the clamp on the BRB 1 is kept open.

REINFUSING SALVAGED BLOOD

⚠ WARNING

- It is recommended that intraoperative and postoperative autologous blood be processed. The washing operation is necessary whenever the recovery is carried out in the presence of heavy haemolysis caused by surgery, or in cases of high coagulopathy or hepatic or renal failure.
- The decision to re-infuse the blood components obtained by this concentration-washing procedure is the responsibility of the attending physician.
- Do not reinfuse the patient's blood from the primary reinfusion bag when it is connected to the cell separator's wash circuit. Reinfusion from the primary reinfusion bag connected to the circuit could lead to air embolism to the patient.
- Do not reinfuse under pressure (i.e., do not use a blood pressure cuff on the secondary reinfusion bag). Reinfusion under pressure could lead to air embolism.
- To minimize the complications of particulate matter infusion, use of an in-line microaggregate filter on the patient reinfusion line is highly recommended. Failure to use an in-line filter may result in reinfusing particles potentially harmful to the patient.
- To reduce the risk of air embolism, remove all air from the primary reinfusion bag before preparing the bag for reinfusion.
- The American Association of Blood Banks recommends the following guidelines for expiration of salvaged blood¹.
 - a. If not transfused immediately, units collected and processed under sterile conditions with a device for intraoperative blood collection that washes with 0.9% saline, USP, shall be stored under one of the following conditions prior to completion of transfusion:
 - At room temperature for up to 4 hours after completion of processing.
 - at 1°C to 6°C for up to 24 hours, provided that storage at 1°C to 6°C is begun within 4 hours after completion of processing
 - b. The transfusion of shed blood collected under postoperative or post-traumatic conditions shall be completed within 6 hours of initiating the collection.

¹ American Association of Blood Banks. Standards for Perioperative Autologous Blood Collection and Administration. 3rd ed. Bethesda, MD: American Association of Blood Banks, 2007.

- Limit the reinfusion of unprocessed blood to the surgical situations of real necessity. The decision to reinfuse unprocessed blood is on the sole responsibility of the physician in charge. If reinfusing unprocessed blood, the use of a microaggregate filter on the patient reinfusion line is mandatory.

04269 WASTE BAG X

⚠ WARNING

The waste bag should be replaced with the equipment in STOP and the bowl completely empty.

SET UP OF WASTE BAG X

1. Hang the waste bag on the two lower hangers located on the right panel of the machine (XTRA). Check that the measuring scale is visible and the drain valve is easily accessed. Prior to using the disposable set, ensure that the drain valve is closed.
2. Ensure that the clamp on the inlet tube of the waste bag is OPEN.
3. If you are not going to perform any preoperative sequestration procedure prior to Intraoperative recovery, connect the waste line from the centrifuge bowl to the waste bag. The connection must be tight.
4. Drainage of the liquid in the bag can be done during the procedure with the equipment in STOP and the bowl completely empty. Turn the valve ring nut anticlockwise and proceed to drain until the level of the liquid reaches the line marked on the bag. Retighten the valve by turning the ring nut clockwise

⚠ CAUTION

- A full waste bag can be emptied by means of the drain valve on the bottom of the bag or replaced aseptically with a new one. Do not completely empty the waste bag until the end of the procedure. If you empty the waste bag during the procedure, leave approximately one liter of fluid in the waste bag to prevent the possibility of vacuum being generated in the waste bag during the Empty cycle. Vacuum in the waste bag may prevent complete emptying of the bowl.
- If clamps are used during bag substitution, make sure to remove the clamps prior to continue blood processing.

⚠ WARNING

- The waste bag should be replaced with the equipment STOPPED and the bowl completely empty.

04270 RESERVOIR OUTLET CLAMP XRES B

SET UP OF THE RESERVOIR OUTLET CLAMP XRES B

1. Using aseptic technique, close slide clamp, remove blue cover and attach to 1/4" outlet port of the XRES B Blood Collection Reservoir.
2. Collect blood into reservoir per hospital protocol.
3. If processing will occur, remove the blue cap from the stepped connector and attach to the Bowl Set. Open slide clamp to initiate flow.
4. Discard after use

04271 RESERVOIR OUTLET "Y" XRES B

SET UP OF THE RESERVOIR OUTLET "Y" XRES B

1. Using aseptic technique, close slide clamps, remove the blue cover and attach 1/4" single female port to the Reservoir Outlet Clamp XRES B.
2. Collect blood into reservoir per hospital protocol.
3. If processing will occur, remove the blue cap from the male connector and attach to the Bowl Set. Open slide clamp to initiate flow.
4. For connection to the ECC circuit, access the 1/4" female port, removing the blue cap. Regulate flow through the port with the slide clamp.
5. Discard after use

MEDICAL DEVICES TO BE USED WITH THE PRODUCTS

To carry out intraoperative and postoperative blood treatments the series of disposable accessories must be used in conjunction with the SORIN GROUP ITALIA medical devices

⚠ CAUTION

The User should observe the warnings and cautions and follow Instruction for Use accompanying the separate device.

RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to bloodborne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned goods authorization (RGA) number is required from Sorin Group USA, Inc. prior to shipping. If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:

SORIN GROUP USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031

LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected. SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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